Magnetostrictive Ultrasonic Inserts

INSTRUCTIONS FOR USE

Intended Purpose:

The Magnetostrictive Ultrasonic Inserts are the working-end devices that remove calculus deposits from teeth and debris from gingiva during ultrasonic dental cleaning and periodontal therapy.

Indications:

The Magnetostrictive Ultrasonic Inserts generate ultrasonic vibrations intended for ultrasonic procedures:

1. General Supra and Subgingival scaling applications.
2. Periodontal debridement for non-surgical treatment of Periodontal diseases.

Contraindications:

The Magnetostrictive Ultrasonic Inserts should not be used in patients with cardiac pacemakers. Devices which show evidence of corrosion should not be used as they may pose an elevated hypersensitivity risk to patients with nickel-allergy. Do not expose to phenols or iodophors. Do not use dry heat sterilization, or heat above 135°C (275°F).

Intended Patient Population:

The intended patient population is adult and pediatric patients under the care of a dental professional.

Description:

The Magnetostrictive Ultrasonic Inserts are the metal working-ends attached to a grip and metal stack that is inserted into the handpiece in any major manufacturers ultrasonic scaling device. For intended use by dental professionals only. The inserts feature a variety of tip designs, water delivery methods, grips, handle colors and frequencies (either 25 kHz or 30 kHz).

General Considerations:

First use instructions:

- After receiving a new Magnetostrictive Ultrasonic Insert, make sure to follow initial cleaning, disinfection and sterilization steps before using for the first time. All non-sterile instruments are to be cleaned, disinfected, and sterilized prior to first use and each subsequent use.

Sterilizer and automatic cleaning equipment best practices:

- Sterilization and cleaning equipment should be maintained and checked per the manufacturer's recommendation.
- Consider the legal provisions valid for your country as well as any hygienic instructions of the doctor's practice or hospital.

Water quality

- The use of completely deionized or distilled water can help avoid most water quality problems.

Instructions For Use:

For proper use of ultrasonic inserts, please review the following instructions:

- Fill the handpiece with water until slightly overflowing.
- Lubricate the O-ring (black or green rubber material) with water.
- Lightly twist the insert down into the handpiece until the insert is firmly in place.

For proper insert selection, please select according to the following:

- Evaluation of Patient health or disease status
- Type and location of deposits or stain to be removed.

Please use the recommended power settings according to tip design per the table below. Always start with the lowest recommended power setting.

Insert Tip	Tip Description	Tip Cross Section	Recommendation	Recommended Power Settings	Frequency	Handle Type and Water Delivery	Handle Color	Grip	Part Codes
#10 Universal	Standard diameter shank extending to a tapered tip.		Light to Heavy Deposits Tenacious Calculus Orthodontic Cement All-purpose insert for supragingival moderate to heavy debris removal. Supragingival use, can be used subgingival if tissue allows.	Low - High	25kHz, 30kHz	Swivel Direct Flow	Lavender with grey cone	Silicone	UI30SD10 UI25SD10
						Swivel Base Flow	Lavender with black cone	Silicone	UI30SS10 UI25SS10
						Streamline Direct Flow	Lavender	Resin	UI30SF10 UI25SF10
						Streamline Base Flow	Lavender	Resin	UI1030KX UI1025K
#100 Thin	Long thin shank extending to a thin tapered tip.		Light to Moderate Deposits Thin shank allows effective subgingival and interproximal access. Ideal for periodontal maintenance and localized periodontal defects.	Low - Med	25kHz, 30kHz	Swivel Direct Flow	Black with grey cone	Silicone	UI30SD100 UI25SD100
						Swivel Base Flow	Black with black cone	Silicone	UI30SS100 UI25SS100
						Streamline Direct Flow	Black	Resin	UI30SF100 UI25SF100
						Streamline Base Flow	Black	Resin	UI30K100SX UI25K100S
#1000 Triple Bend	Triple bend shank extending to a tapered tip with beveled edges.		Moderate to Heavy Deposits Tenacious Deposits Orthodontic Cement Triple bend shank allows increased access around line angles and interproximally. Beveled tip concentrated energy at corners allowing for enhanced calculus and stain removal. Supragingival use, subgingival if tissue allows.	Low - High	25kHz, 30kHz	Swivel Direct Flow	Orange with grey cone	Silicone	UI30SD1000 UI25SD1000
						Swivel Base Flow	Orange with black cone	Silicone	UI30SS1000 UI25SS1000
						Streamline Direct Flow	Orange	Resin	UI30SF1000 UI25SF1000
						Streamline Base Flow	Orange	Resin	UI100030KX UI100025K

Insert Tip	Tip Description	Tip Cross Section	Recommendation	Recommended Power Settings	Frequency	Handle Type and Water Delivery	Handle Color	Grip	Part Codes
XT Triple Bend	Thin triple bend shank extending to a thin beveled edge tip 24% thinner than #1000 insert		Light to Moderate Deposits Thin triple bend shank allows increased access around line angles and interproximally. Beveled tip concentrated energy at corners allowing for enhanced calculus and stain removal. For both Supra and subgingival surfaces.	Low - Med	30kHz	Swivel Direct Flow	Purple	Silicone	UI30SDXTTB
						Streamline Direct Flow	Purple	Resin	UI30SFXTTB
XT	Long, ultra thin shank tapers to a very thin tip.		Light to Moderate Deposits Tenacious Subgingival Deposits Site specific uses including narrow, deep pockets, furcations, developmental concavities and tight areas between misaligned teeth.	Low - Med	25kHz, 30kHz	Swivel Base Flow	Purple with black cone	Silicone	UI30SSXT UI25SSXT
						Streamline Direct Flow	Purple	Resin	UI30SFXT
After Five, Left	Extended thin shank narrows to a tapered, probe sized tip.		Light to Moderate Deposits Allows improved access to posterior proximal surfaces. Can be used in furcations. Available alone or as a set with right and straight inserts.	Low - Med	25kHz, 30kHz	Swivel Direct Flow	Green	Silicone	UI30SDL UI25SDL
						Streamline Direct Flow	Green	Resin	UI30SFL UI25SFL
						After Five Plus External Water Tube	Green	Resin	UI30KSF10L UI25KSF10L
						After Five External Water Tube	Metal	Metal	UI30KSL10L UI25KSL10L
After Five, Right	Extended thin shank narrows to a tapered, probe sized tip.		Light to Moderate Deposits improved access to posterior proximal surfaces. Can be used in furcations. Available alone or as a set with left and straight inserts.	Low - Med	25kHz, 30kHz	Swivel Direct Flow	Red	Silicone	UI30SDR UI25SDR
						Streamline Direct Flow	Red	Resin	UI30SFR UI25SFR
						After Five Plus External Water Tube	Red	Resin	UI30KSF10R UI25KSF10R
						After Five External Water Tube	Metal	Metal	UI30KSL10R UI25KSL10R
After Five, Straight	Extended thin shank narrows to a tapered, probe sized tip. 40% thinner than original prophyl design.		Light to Moderate Deposits Used to remove calculus from periodontal pockets with improved adaptation for deeper pockets. Available alone or as a set with left and right inserts	Low - Med	25kHz, 30kHz	Swivel Direct Flow	Blue	Silicone	UI30SDS UI25SDS
						Streamline Direct Flow	Blue	Resin	UI30SFS UI25SFS
						After Five Plus External Water Tube	Blue	Resin	UI30KSF10S UI25KSF10S
						After Five External Water Tube	Metal	Metal	UI30KSL10S UI30KSL10S
#3 Beaver-tail	Thick, broad, flat tip.		Heavy supragingival deposits Ideal for anterior lingual surfaces, tenacious stain, orthodontic cement	Med - High	25kHz, 30kHz	Streamline Base Flow	Yellow	Resin	UI330KX UI325K
P10 & P10P	Standard diameter shank extending to a thin tapered tip. Original prophyl style inserts with external water tube.		Light to Heavy deposits. Orthodontic cement Supragingival and subgingival.	Low- High	25kHz, 30kHz	Original Prophyl External Water Tube	Metal	Metal	UI30KP10 UI25KP10
						Prophyl Plus External Water Tube	Green	Resin	UI30KP10P UI25KP10P
Power Plus Standard Conical	Standard diameter shank extending to an extended length tapered tip.		Moderate to heavy deposits. Supragingival, wide, deep pockets.	Low - High	30kHz	Swivel Base Flow	Teal with black cone	Resin	UI30PPSC

Insert Tip	Tip Description	Tip Cross Section	Recommendation	Recommended Power Settings	Frequency	Handle Type and Water Delivery	Handle Color	Grip	Part Codes
Power Plus Thin Conical	Thin diameter shank extending to an extended length tapered tip.		Light to moderate deposits. Narrow, deep vertical pockets. Areas mesial to retromolar pad, crowded areas, general supra use.	Low – Med	30kHz	Swivel Base Flow	Teal with black cone	Resin	UI30PPTC
Power Plus Standard Beveled	Standard diameter shank extending to an extended length tapered tip.		Moderate to heavy deposits. Tenacious deposits. Supragingival, wide, deep pockets.	Low - High	30kHz	Swivel Base Flow	Teal with black cone	Resin	UI30PPSB
Power Plus Thin Beveled	Thin diameter shank extending to an extended length tapered tip.		Light to moderate deposits. Tenacious deposits. Narrow, deep vertical pockets, crowded areas.	Low – Med	30kHz	Swivel Base Flow	Teal with black cone	Resin	UI30PPTB

The metal tips of the ultrasonic insert will wear with each additional use. Tips should be regularly measured since tip wear affects insert performance. The insert wear guide can be used to determine efficiency of the insert. A loss of 1 mm (blue line) represents a 25% loss of efficiency. A loss of 2 mm (red line) represents a 50% loss of efficiency, and the insert should be replaced.

Warnings:

The Magnetostrictive Ultrasonic Insert tips (i) should be replaced when "2mm wear occurs"; (ii) finger stick/cut, scrapped finger on insert, cut hand/finger, or needle stick/puncture may occur with potential for the transfer of bloodborne pathogens during tip removal from handpiece if care is not exercised; (iii) 'material fragmentation' (e.g. broken tip) may occur if used improperly; (iv) accidental swallowing of inserts/tips/tip fragments may occur; (v) device fracture/break may occur if used improperly; (vi) foreign body in patient may occur (e.g. tip swallowed); (vii) device embedded in tissue or plaque may occur; (viii) perforation of tooth may occur; (ix) patient discomfort (e.g pain, burn or bleeding) may occur; (x) exposure to blood/bodily fluids; (xi) patient infection/Sepsis may occur; (xii) Unintentional electrical shock from the generator (unit) may occur; (xiii) leaking may occur; and (xiv) loose housing may occur.

Precautions:

- Use aerosol mitigation protocols (pre-treatment with antibacterial rinses, high-volume evacuator, etc.) and proper PPE (surgical mask covering nose and mouth, gowns/jackets in work area, medical grade single-use gloves, etc.) to protect from infection and burns.
- Wash hands before and after treating each patient, before putting on gloves, and immediately after removing gloves.
- Do not disassemble or alter the insert.
- Check insert vibration outside the oral cavity before use.
- Do not exceed recommended power ranges for the inserts. Excess power may damage tooth structure and/or result in damage to or excess wear of the insert.
- Check the insert before every use. Inserts that have been bent, altered, worn or compromised in any way should be removed from service.
- Do not use inserts on metal or ceramic restorations, unless inserts are intended specifically for this purpose.
- To eliminate risk of damage to the soft and hard tissues keep the insert tip moving at all times when in contact with the tooth, gingiva, mucosa and intra-oral tissues.
- Retract lips, cheeks and tongue to prevent contact with the working end of the insert
- Make sure the green or black O-ring within the insert is in-place prior to use, and lubricate with water prior to inserting into the handpiece.
- Do not apply the insert tip directly to the tooth surface.

Manual Cleaning & Disinfection

To clean, rinse inserts thoroughly or fully immerse in a mild ultrasonic cleaning solution. Using a soft brush, brush inserts to remove residues from the surface. Inserts may be placed in an ultrasonic cleaning unit for 7-10 minutes or a minimum of 16 minutes if using an IMS™ Cassette.

Thoroughly rinse with tap water and dry completely prior to sterilization. Do not apply excessive force during cleaning and reprocessing.

Automated Cleaning & Disinfection

Load washer-disinfector as validated for the intended cycle, being careful not to overload. Start the cycle. Ensure instruments are dry prior to handling. Remove instruments following cycle completion.

The fundamental suitability of the instruments for effective automated cleaning and disinfection was demonstrated by an independent accredited test laboratory under the following conditions:

Washer-Disinfector	Miele Professional G 7836 CD
Racks	Mobile injector unit (Miele) E429, Four-level rack (Miele) E 493
Cleaning Cycle	2 minute pre-cleaning with cold tap water, Draining 5 minute cleaning with 55 °C cleaning solution, Draining 3 minute rinsing with cold deionized water, Draining 2 minute rinsing with cold deionized water, Draining
Cleaning Solution	0.5 % cleaning solution neodisher® Mediclean Dental (Chemische Fabrik Dr. Weigert, Hamburg)
Validation Report	Project Number: 00418-1 23919 Examination of an Automated Cleaning Process using quantitative Detection of Protein and Hemoglobin and the Radionuclide Method

Sterilization

Sterilize inserts in an IMS™ cassette or paper pouch in a steam autoclave. Steam sterilize for at least 4 minutes at 270°F/132°C or 30 minutes at 250°F/121°C. Do not heat above 275°F/135°C.

30 minutes minimum of drying time is recommended after steam sterilization. Store sterilized inserts in a dry, dust free location. Sterilization can only be maintained if the inserts remain packaged or wrapped.

Precautions:

Do not use chemical or surface disinfectants on ultrasonic inserts as deterioration of the material may occur. Avoid contact with corrosive chemicals such as chlorides, sulfates and surface disinfectants. Do not expose the insert to solutions containing phenols or iodophors. Chemical vapor sterilization, rapid heat transfer sterilization, and cold liquid/chemical sterilants are also not recommended and may void the warranty.

The fundamental suitability of the instruments for an effective sterilization was demonstrated by an independent accredited test laboratory under the following conditions:

Sterilization Method	Pre-vacuum Mode
Sterilizer	W & H Lisa MB 17 Steam Sterilizer
Sterilization Temperature	134 °C (273°F)
Pre-Vacuum Phases	3
Holding (full cycle)	4 minutes
Drying Time	30 minutes
Validation Report	Project Numbers: 25517-1; 25517-2; 24219 Validation of a Sterilization Process using Steam Sterilization in Pre-vacuum Mode Method MD 4.0: Sterilization validation of medical products with moist heat Project Numbers: 10918-1; 10918-2 Determination of Residual Moisture after Sterilization using Steam Sterilization in Pre-vacuum Mode

Compatible Magneto Units:


Manufacturer	Unit
DENTSPLY	Bobcat Pro 25k
DENTSPLY	Cavitron Select SPS
DENTSPLY	Cavitron Plus
DENTSPLY	Cavitron 300 Series
PRIVATE LABEL	Schein Acclean
COLTENE	Biosonic US100R
PARKELL	Turbo Sensor

Disposal:


Handle and dispose of broken tips in appropriate waste receptacle per facility protocols. Dispose inserts per local regulation.

Ultrasonic Insert Wear & Power Level Guide

1. Position point A to point B as shown.
2. Rotate insert until tip is flat against the card.
3. Evaluate tip condition, compared to gray/red bars.

Optimum Efficiency
 25% Efficiency Loss (Reorder)
 50% Efficiency Loss (Discard)

In case of a serious incident that has occurred in relation to Hu-Friedy instruments within the European Union, please report by either calling Hu-Friedy for assistance Tel 00800 4837 4339 or send an e-mail to info@hu-friedy.eu. In addition, please inform your national competent authority.

Hu-Friedy Mfg. Co., LLC
 2501 Sanders Rd.
 Northbrook, IL 60062 | USA
 1-800-Hu-Friedy | HuFriedyGroup.com

EU REP Hu-Friedy Europe LLC & Co. KG
 Bogenstrasse 6 - 8
 78576 Emmingen - Liptingen
 Germany

CH REP MDSS CH GmbH
 Laurenzenvorstadt 61
 5000 Aarau
 Switzerland

UK REP MDSS-UK RP Ltd.
 Parkway House, Palatine Rd,
 Northenden, Wythenshawe,
 Manchester M22 4DB
 United Kingdom

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
 ©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved.

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©2026 Hu-Friedy Mfg. Co., LLC. Alla rättigheter förbehållna.